

# Revêtement multicouche hautement hermétique pour implants miniaturisés

*Aperçu de la technologie*

Hicham Damsir

Medi'nov connection Lyon Avril 2019



1. Comelec in key figures & activities
2. Challenges for electronic implantable devices
3. Technology overview
4. Case study 1 : Implanted magnet used in miniaturized pumping system
5. Case study 2 : Wireless implanted sensor
6. Conclusions

# Comelec in key figures & activities

**40 years** experience  
in Parylene



Member of  
watchmaking group  
**~100 people**



**2000** run  
per year



**60%** to the export

Equipment Sales over  
**3 continents**

- Parylene and customized **coating service**
  - Parylene N,C, D, F-VT4, F-AF4
  - Coating based PECVD and ALD ( $\text{Al}_2\text{O}_3$ ,  $\text{SiO}_2$ ,  $\text{TiO}_2$ )
- Masking/Demasking service
- Parylene etching technologies
- **Equipment manufacturing** & technological transfer
- Distribution of raw material and consumables
- Consulting on specific demand (clean room integration, new process, material,...)



# Comelec in key figures & activities



# EU R&D Program

- ✓ Innosuisse, Swiss funded projects
- ✓ European projects



2007

2010

2015

2018

2021

**Deep brain stimulator:** Electrode protection

**Therapeutic Bioelectrode:** Coating of PCBA

**Valve:** Coating of sensitive magnetic alloy

**Tips:** Biocompatible coating


**Bone implant sensors:** Protective coating

**Glucose sensors:** Coating of PCBA

**Liquid pumping systems:** Coating of PCBA driver

**Earing aids:** Coating of PCBA

**RIFD chip:** Glass tag coating



**Artificial heart:** Protective layer from blood contact

**Heart pump:** Protective layer and encapsulation

**Pacemaker:** coating of the PCB


**Silicone prosthesis:** Anti-friction coating

**Rubber piston:** anti-sticking

- Miniaturization leads to high density of components per surface unit and a **risk of current leakage**.
- Development of MEMS and microfluidic systems enable a fast and precise diagnostic. However their **sensitivity to particles** and physiological fluids decrease dramatically their effectiveness.
- For a fast connectivity, the component used (eg : RF antenna) **thin packaging** is required.
- In advanced electronic implants, several material **need protection against leaching**.
- For long term implant (More than 1 year), the **resistance to highly conductive and corrosive physiological fluids** becomes a crucial issue.
- Bacterial colonization is also a topic of great concern

**The encapsulation layer is key component to overcome these challenges!**






**Parylene C 25  $\mu\text{m}$  :**  
WVTR => ~ 7  $\text{g}/\text{m}^2/\text{day}$   
O2TR=> ~ 120  $\text{cm}^3/\text{m}^2/\text{day}$


Targeted performance for long term AIMD

## Why combining Parylene with oxide material? The synergetic effect...

- Parylene C barrier properties are limited
- The oxide material lead to weak point that let molecules diffuse
- The combination lead to full encapsulation of the oxide material and stress dissipation.
- The single sandwich can be repeated many times if necessary



Quantitative result : OTR, WTR



## Context :

- Magnet are made of Neodymium Iron Bore Alloy which is extremely sensitive to conductive liquids
- The product is in direct contact with blood and other potential physiological fluids
- **Magnet need to resist over 50 days in salt solution at room temperature without any sign of corrosion**



*Example of magnet*

## Test results :

- A multilayer structure has been proposed based on ALD Aluminum oxide and Parylene C.
- The multilayer was compared with the current solution based on Parylene C 15  $\mu\text{m}$  and a single ALD layer



- ✓ ALD single layer
- ✓ Corrosion after 2 days



- ✓ Parylene C single layer
- ✓ No corrosion after 50 days
- ✓ Some microscopic spot could be observed




- ✓ Parylene C ALD multilayer
- ✓ No corrosion after 50 days

## Case study 2 : Wireless implanted glucose sensor


### Context :

- The traditional **metallic cap are not allowed** due the presence of antenna and the need to access directly to the protein.
- A **thin encapsulation layer is needed** to do not disturb the signal and the RF communication.
- The implant need to stay operational **more than 1 year**
- A current solution based on Parylene 20  $\mu\text{m}$  is working for up to 3 months.



**Test results :**

- To extend the device shelf life, a multilayer coating based on Al<sub>2</sub>O<sub>3</sub>/TiO<sub>2</sub> + Parylene C 22 μm is proposed
- According to simulation, an extension of the shelf life by three years is expected.



Structure	Time where a defect is detected After aging 55°C in PBS	Extrapolation @ 37°C in vivo
Parylene C 22 μm	3 weeks	3 months
Parylene C/Al <sub>2</sub> O <sub>3</sub> /TiO <sub>2</sub>	12 weeks* (no defect detected, test under progress)	12 months*
PVDTiO <sub>2</sub> / ALD Al <sub>2</sub> O <sub>3</sub> /TiO <sub>2</sub> /Parylene C	6 weeks	6 months

- This **multilayer high barrier technology** open the door for new range of application where thickness/permeability ratio is critical **in AIMD** especially when sensors are involved
- It keeps the same advantages of the well known Parylene conformal coating at a **thinner layer**
- Performance can be tuned by increasing the **number of layers**
- **Biocompatibility** can be easily demonstrated thanks to Parylene top coating
- **Cost effective process** based on well-established technologies
- **Comelec SA Parylene Coating is ready for new cooperation projects**







Hicham Damsir  
Mat Sc. Eng  
Business Manager  
Headquarter

Tel +41329240004

[h.damsir@comelec.ch](mailto:h.damsir@comelec.ch)

**THANK YOU FOR YOUR ATTENTION!**